

SHINVA



CIMA[®]
Industries Inc.
Pharmaceutical Equipment

PRESSURE STEAM STERILIZER ⚙️



*Integral solutions
in medical equipment*

PRESSURE STEAM STERILIZER

USAGE

It is applicable for the sterilization of unpacked solid instrument in department of stomatology, ophthalmology, operation room and CSSD.

- Liquid-crystal display of temperature, pressure, time, running status and failure warning.
- Equipped with professional and systematized sterilizing programs fabric, instrument, rubber and liquid.
- The positive pressure air-discharging system will effectively discharge the cold air from chamber to improve steam saturability.
- Build-in efficient steam generator, microcomputer temperature control, time and energy saving.
- Equipped with post-drying function which makes the residual humidity of sterilized articles less than 3%.
- Equipped with printer to print date, time and process parameters.
- Over-heat auto protection device, double-way door safety lock device, overpressure relief valve, short-circuit protection and overpressure protection.
- Perfect four-level authority management system to avoid misoperation.



XG1 .CD-135M
XG1 .CD-185M
XG1 .CD-300M



XG1 .CD-60M
XG1 .CD-80M

Main technical parameters

- Design pressure: 0.28MPa
- Maximum working pressure: 0.23MPa
- Sterilizing temperature: 105-136 e
- Chamber material: 06Cr19Ni10 stainless steel (304)
- Accuracy of temperature display: 0.1°C
- Accuracy of pressure display: 1 Kpa

Order Code	Volume (L)	Chamber Dimension (φxL)	Overall Dimension (LxWxD) mm	Net Weight	Voltage	Power Frequency	Power
XGI .CD-60M	60L	380x510 mm	1000x697x1600 mm	230Kg	380V	50Hz	10Kw
XGI .CD-80M	80L	380x710 mm	1000x697x1600 mm	260Kg	380V	50Hz	10Kw
XGI .CD-135M	135L	500x710 mm	1110x750x1742 mm	400Kg	380V	50Hz	10Kw
XGI .CD-185M	185L	500x950 mm	1350x750x1742 mm	460Kg	380V	50Hz	10Kw
XGI .CD-300M	300L	632x1000 mm	1400x890x1780 mm	750Kg	380V	50Hz	16Kw



World Trade Center. 1er Piso, Área Comercial
Calle 53 Marbella. Apartado 0832-00155 WTC
Panamá, República de Panamá
Tel. +507 205 1915 Fax. +507 205 1802

Address: No. 2009 Xinglu Ave. Shinva Science Park, Zhoucun, Zibo,
Shandong Province, P. R. China
Postcode:255300
Sales telephone: 0086-533-3581273
Email: xhmed@vip.163.com
Website: www.shinva.com